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1. Hoffman M, Monroe DM, 3rd. A cell-based model of hemostasis. *Thromb Haemost* 2001;85:958–65.
2. Hartert H. Thrombelastography: physical and physiological aspects. In: Copley AL, Stainsby G, eds. *Flow properties of blood and other biological systems*: Pergamon Press, 1960:186–96.
3. Nair SC, Dargaud Y, Chitlur M et al. Tests of global haemostasis and their applications in bleeding disorders. *Haemophilia* 2010;16 Suppl 5:85–92.
4. Tanaka KA, Szlam F, Sun HY et al. Thrombin generation assay and viscoelastic coagulation monitors demonstrate differences in the mode of thrombin inhibition between unfractionated heparin and bivalirudin. *Anesth Analg* 2007;105:933–9, table of contents.
5. Cotton BA, Faz G, Hatch QM et al. Rapid thrombelastography delivers real-time results that predict transfusion within 1 hour of admission. *J Trauma* 2011;71:407–14.

6. Kashuk JL, Moore EE, Sabel A et al. Rapid thrombelastography (r-TEG) identifies hypercoagulability and predicts thromboembolic events in surgical patients. *Surgery* 2009;146:764–72.
7. Riskin DJ, Tsai TC, Riskin L et al. Massive transfusion protocols: the role of aggressive resuscitation versus product ratio in mortality reduction. *J Am Coll Surg* 2009;209:198–205.
8. Sambu N, Curzen N. Monitoring the effectiveness of antiplatelet therapy: opportunities and limitations. *Br J Clin Pharmacol* 2011;72:683–96.
9. Kang Y. Thromboelastography in liver transplantation. *Semin Thromb Hemost* 1995;21 Suppl 4:34–44.
10. Kozek-Langenecker SA, Afshari A, Albaladejo P et al. Management of severe perioperative bleeding: Guidelines from the European Society of Anaesthesiology. *Eur J Anaesthesiol* 2013;30:270–382.
11. Wang SC, Shieh JF, Chang KY et al. Thromboelastography-guided transfusion decreases intraoperative blood transfusion during orthotopic liver transplantation: randomized clinical trial. *Transplant Proc* 2010;42:2590–3.
12. National Institute for Health and Care Excellence. Detecting, managing and monitoring haemostasis: viscoelastometric point-of-care testing (ROTEM, TEG and Sonoclot systems). <https://www.nice.org.uk/guidance/dg13/resources/nice-diagnostic-support-for-viscoelastometric-pointofcare-testing-rotem-teg-and-sonoclot-systems-99757>. 2014; Accessed April 2015.

13. American Society of Anesthesiologists Task Force on Perioperative Blood Management. Practice guidelines for perioperative blood management: an updated report by the American Society of Anesthesiologists Task Force on Perioperative Blood Management*. *Anesthesiology* 2015;122:241–75.
14. Shore-Lesserson L, Manspeizer HE, DePerio M et al. Thromboelastography-guided transfusion algorithm reduces transfusions in complex cardiac surgery. *Anesth Analg* 1999;88:312–9.
15. Ak K, Isbir CS, Tetik S et al. Thromboelastography-based transfusion algorithm reduces blood product use after elective CABG: a prospective randomized study. *J Card Surg* 2009;24:404–10.
16. Nuttall GA, Oliver WC, Santrach PJ et al. Efficacy of a simple intraoperative transfusion algorithm for nonerythrocyte component utilization after cardiopulmonary bypass. *Anesthesiology* 2001;94:773–81.
17. Agarwal S, Johnson RI, Shaw M. Preoperative point-of-care platelet function testing in cardiac surgery. *J Cardiothorac Vasc Anesth* 2015;29:333–41.

18. Hendriks HG, van der Maaten JM, de Wolf J et al. An effective treatment of severe intractable bleeding after valve repair by one single dose of activated recombinant factor VII. *Anesth Analg* 2001;93:287–9.
19. Alfirevic A, Duncan A, You J et al. Recombinant factor VII is associated with worse survival in complex cardiac surgical patients. *Ann Thorac Surg* 2014;98:618–24.
20. Hendriks HG, Meijer K, de Wolf JT et al. Effects of recombinant activated factor VII on coagulation measured by thromboelastography in liver transplantation. *Blood Coagul Fibrinolysis* 2002;13:309–13.
21. Kawaguchi C, Takahashi Y, Hanesaka Y et al. The in vitro analysis of the coagulation mechanism of activated factor VII using thrombelastogram. *Thromb Haemost* 2002;88:768–72.
22. Qi X, Zhao Y, Li K et al. Evaluating and monitoring the efficacy of recombinant activated factor VIIa in patients with haemophilia and inhibitors. *Blood Coagul Fibrinolysis* 2014;25:754–60.
23. Hunt H, Stanworth S, Curry N et al. Thromboelastography (TEG) and rotational thromboelastometry (ROTEM) for trauma-induced coagulopathy in adult trauma patients with bleeding. *Cochrane Database Syst Rev* 2015;2:CD010438.

24. Spahn DR, Bouillon B, Cerny V et al. Management of bleeding and coagulopathy following major trauma: an updated European guideline. *Crit Care* 2013;17:R76.
25. Thomas D, Wee M, Clyburn P et al. Blood transfusion and the anaesthetist: management of massive haemorrhage. *Anaesthesia* 2010;65:1153–61.
26. American College of Surgeons Trauma Quality Improvement Program. Massive transfusion in trauma guidelines.
<https://www.facs.org/~media/files/quality%20programs/trauma/tqip/massive%20transfusion%20in%20trauma%20guidelinesashx>; Accessed April 2015.
27. Johansson PI. Coagulation monitoring of the bleeding traumatized patient. *Curr Opin Anaesthesiol* 2012;25:235–41.
28. Johansson PI, Sorensen AM, Larsen CF et al. Low hemorrhage-related mortality in trauma patients in a Level I trauma center employing transfusion packages and early thromboelastography-directed hemostatic resuscitation with plasma and platelets. *Transfusion* 2013;53:3088–99.
29. Johansson PI, Stensballe J, Oliveri R et al. How I treat patients with massive hemorrhage. *Blood* 2014;124:3052–8.

30. Kashuk JL, Moore EE, Wohlauer M et al. Initial experiences with point-of-care rapid thrombelastography for management of life-threatening postinjury coagulopathy. *Transfusion* 2012;52:23–33.
31. Kane I, Ong A, Orozco FR et al. Thromboelastography predictive of death in trauma patients. *Orthop Surg* 2015;7:26–30.
32. Johansson PI. The blood bank: from provider to partner in treatment of massively bleeding patients. *Transfusion* 2007;47:176S–81S.
33. Gonzalez E, Moore EE, Moore HB et al. Trauma-induced coagulopathy: An institution's 35 year perspective on practice and research. *Scand J Surg* 2014;103:89–103.
34. Holcomb JB, Minei KM, Scerbo ML et al. Admission rapid thrombelastography can replace conventional coagulation tests in the emergency department: experience with 1974 consecutive trauma patients. *Ann Surg* 2012;256:476–86.
35. Kashuk JL, Moore EE, Sawyer M et al. Postinjury coagulopathy management: goal directed resuscitation via POC thrombelastography. *Ann Surg* 2010;251:604–14.
36. Sawyer MM, Myers G, Humphrey J et al. Trauma and thrombelastography: how changes in the understanding of coagulopathy, testing, and hospital systems have changed one group's practice. *Semin Cardiothorac Vasc Anesth* 2012;16:142–52.

37. Yin J, Zhao Z, Li Y et al. Goal-directed transfusion protocol via thrombelastography in patients with abdominal trauma: a retrospective study. *World J Emerg Surg* 2014;9:28.
38. Tapia NM, Chang A, Norman M et al. TEG-guided resuscitation is superior to standardized MTP resuscitation in massively transfused penetrating trauma patients. *J Trauma Acute Care Surg* 2013;74:378-85.
39. Gurbel PA, Bliden KP, Guyer K et al. Platelet reactivity in patients and recurrent events post-stenting: results of the PREPARE POST-STENTING Study. *J Am Coll Cardiol* 2005;46:1820-6.
40. McCrath DJ, Cerboni E, Frumento RJ et al. Thromboelastography maximum amplitude predicts postoperative thrombotic complications including myocardial infarction. *Anesth Analg* 2005;100:1576-83.
41. Schreiber MA, Differding J, Thorborg P et al. Hypercoagulability is most prevalent early after injury and in female patients. *J Trauma* 2005;58:475-80.
42. Cotton BA, Minei KM, Radwan ZA et al. Admission rapid thrombelastography predicts development of pulmonary embolism in trauma patients. *J Trauma Acute Care Surg* 2012;72:1470-5.

43. Mahla E, Suarez TA, Bliden KP et al. Platelet function measurement–based strategy to reduce bleeding and waiting time in clopidogrel–treated patients undergoing coronary artery bypass graft surgery: the timing based on platelet function strategy to reduce clopidogrel–associated bleeding related to CABG (TARGET–CABG) study. *Circ Cardiovasc Interv* 2012;5:261–9.
44. Kasivisvanathan R, Abbassi–Ghadi N, Kumar S et al. Risk of bleeding and adverse outcomes predicted by thromboelastography platelet mapping in patients taking clopidogrel within 7 days of non–cardiac surgery. *Br J Surg* 2014;101:1383–90.
45. Hobson AR, Agarwala RA, Swallow RA et al. Thrombelastography: current clinical applications and its potential role in interventional cardiology. *Platelets* 2006;17:509–18.
46. Ferraris VA, Ferraris SP, Moliterno DJ et al. The Society of Thoracic Surgeons practice guideline series: aspirin and other antiplatelet agents during operative coronary revascularization (executive summary). *Ann Thorac Surg* 2005;79:1454–61

47. Lipets EN, Ataulakhanov FI. Global assays of hemostasis in the diagnostics of hypercoagulation and evaluation of thrombosis risk. *Thromb J* 2015;13:4.
48. Krzanicki D, Sugavanam A, Mallett S. Intraoperative hypercoagulability during liver transplantation as demonstrated by thromboelastography. *Liver Transpl* 2013;19:852–61.

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Copyright 2015、Schattauer の許諾による]